



Recombinant Human WD repeat-containing protein 74 (WDR74)

Product Code	CSB-BP747557HU
Abbreviation	WDR74
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6RFH5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAARWNH VVVGTTETGIL KGVNLQRKQA ANFTAGGQPR REEAVSALCW GTGGETQMLV GCADRTVKHF STEDGIFQGQ RHCPGGEGMF RGLAQADGTL ITCVDSGILR VWHDKDKDTS SDPLLELRVG PGVCRMRQDP AHPHVVATGG KENALKIWDL QGSEEPVFRA KNVRNDWLDL RVPIWDQDIQ FLPGSQKLV CTGYHQVRVY DPASPQRRPV LETTYGEYPL TAMTLTPGGN SVIVGNTHGQ LAEIDLRQGR LLGCLKGLAG SVRGLQCHPS KPLLASCGLD RVLRIHRIQN PRGLEHKVYL KSQLNCLLLS GRDNWEDEPQ EPQEPNKVPL EDTETDELWA SLEAAKRKL SGLEQPQ GAL QTRRRKKKRP GSTSP
Source	Baculovirus
Target Names	WDR74
Protein Names	Recommended name: WD repeat-containing protein 74 Alternative name(s): NOP seven-associated protein 1
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.