



# Recombinant Human WD repeat-containing protein 74 (WDR74)

<b>Product Code</b>	CSB-MP747557HU
<b>Abbreviation</b>	WDR74
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6RFH5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAAAARWNH VWVGTETGIL KGVNLQRKQA ANFTAGGQPR REEAVSALCW GTGGETQMLV GCADRTVKHF STEDGIFQGQ RHCPGGEGMF RGLAQADGTL ITCVDSGILR VWHDKDKDTS SDPLLELRVG PGVCRMRQDP AHPHVVATGG KENALKIWDL QGSEEPVFRA KNVRNDWLDL RVPIWDQDIQ FLPGSQKLV CTGYHQVRVY DPASPQRRPV LETTYGEYPL TAMTLTPGGN SVIVGNTHGQ LAEIDLRQGR LLGCLKGLAG SVRGLQCHPS KPLLASCGLD RVLRIHRIQN PRGLEHKVYL KSQLNCLLLS GRDNWEDEPQ EPQEPNKVPL EDTETDELWA SLEAAKRKL SGLEQPQ GAL QTRRRKKKRP GSTSP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WDR74
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 74 Alternative name(s): NOP seven-associated protein 1
<b>Expression Region</b>	1-385
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.