



# Recombinant Rat DnaJ homolog subfamily B member 11 (Dnajb11)

<b>Product Code</b>	CSB-MP747601RA
<b>Abbreviation</b>	Dnajb11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6TUG0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GRDFYKIL GVPRSASIKD IKKAYRKLAL QLHPDRNPDD PQAQEKFQDL GAAYEVLSDS EKRKQYDTYG EEGLKDGHQ S SHGDIFSHFF GDFGFMFGGA PRQQDRNIPR GSDIIVDLEV TLEEVYAGNF VEVVRNKPVA RQAPGKRKCN CRQEMRTTQL GPGRFQMTQE VVCDECPNVK LVNEERTLEV EIEPGVRDGM EYFIGEGEP HVDGEPGDLR FRIKVVKHRI FERRGDDLYT NVTVSLVEAL VGFEMDITHL DGHKVVHISR D KTRPGAKLW KKGEGLPNFD NNNIKGSLII TFDVFDPKEQ LTEEAKEGIK QLLKQGPVQK VYNGLQGY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Dnajb11
<b>Protein Names</b>	Recommended name: DnaJ homolog subfamily B member 11 Alternative name(s): ER-associated DNAJ ER-associated Hsp40 co-chaperone ER-associated dnaJ protein 3 Short name= ERdj3 Short name= ERj3p Liver regeneration-related pro
<b>Expression Region</b>	23-358
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.