



Recombinant Human BRCA1-A complex subunit Abraxas (FAM175A)

Product Code	CSB-EP747632HU-B
Abbreviation	FAM175A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UWZ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGESTSAVL SGFVLGALAF QHLNTDSDTE GFLLGVEVKGE AKNSITDSQM DDVEVVYTID IQKYIPCYQL FSFYNSSGEV NEQALKKILS NVKKNVVGWY KFRRHSDQIM TFRERLLHKN LQEHFSNQDL VFLLLTSPSII TESCSTHRLE HSLYKPQKGL FHRVPLVVAN LGMSEQLGYK TVSGSCMSTG FSRVQTHSS KFFEEDGSLK EVHKINEMYA SLQEELKSIC KKVEDSEQAV DKLVKDVNRL KREIEKRRGA QIQAAAREKNI QKDPQENIFL CQALRTFFPN SEFLHSCVMS LKNRHVSKSS CNYNHHLDVV DNLTLMVEHT DIPEASPAST PQIIKHKALD LDDRWQFKRS RLLDTQDKRS KADTGSSNQD KASKMSSPET DEEIEKMKGF GEYSRSPTF
Source	E.coli
Target Names	ABRAXAS1
Protein Names	Recommended name: BRCA1-A complex subunit Abraxas Alternative name(s): Coiled-coil domain-containing protein 98 Protein FAM175A
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.