



Recombinant Human WD repeat-containing protein 82 (WDR82)

Product Code	CSB-BP747639HU
Abbreviation	WDR82
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UXN9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKLTDSVLRS FRVAKVFREN SDKINCFDFS PNGETVISSS DDDSIVLYDC QEGKPKRTLY SKKYGVDLIR YTHAANTVVY SSKIDDTIR YLSLHDNKYI RYFPGHSKRV VALSMSPVDD TFISGSLDKT IRLWDLRSPN CQGLMHLQ GK PVCSEFDPEGL IFAAGVNSEM VKLYDLRSFD KGP FATFKMQ YDRTCEWTGL KFSNDGKLIL ISTNGSFIRL IDAFKGVVMH TFGGYANSKA VTLEASFTPD SQFIMIGSED GKIHVWNGES GIKVAVLGK HTGPITCLQF NPKFMTFASA CSNMAFWLPT IDD
Source	Baculovirus
Target Names	WDR82
Protein Names	Recommended name: WD repeat-containing protein 82 Alternative name(s): Swd2
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.