



Recombinant *Gossypium hirsutum* Tubulin beta-7 chain

Product Code	CSB-MP747650GHB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6VAF5
Product Type	Recombinant Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Purity	>85% (SDS-PAGE)
Sequence	MREILHVQAG QCGNQIGGKF WEVVSDEHGI DPKGNYVGTS RVQLERVNVY YNEASGGRYV PRAVLMDLEP GTMDSLRTGP HGQLFRPDNF IFGQNGAGNN WAKGHYTEGA ELIDSVLDVV RKEAENCACL QGFQICHSLG GGTGSGMGTL LISKIKEEYP DRMMLTFSVF PSPKVSDTVV EPYNATLSVH QLVENGDECM VLDNEALYDI CFRTLKLTNP SFGDLNRLIS TTMSGATCCL RFPGQLNSDL RKLAVNLIPF PRLHFFMVG F APLTSSSSSQQ YRALTPELT QQMWDARNMM CAADPRHGRY LTASAMFRGK MSTKEVDEQM INVQNKNSSY FVEWIPNNVK SSVCDIPPTG LTMSSTFMGN STSIQEMFRR VSEQFTVMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVS EYQQYQDAVA DDNDEYDE AMEN
Source	Mammalian cell
Protein Names	Recommended name: Tubulin beta-7 chain Alternative name(s): Beta-7-tubulin
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.