



Recombinant *Oryza sativa* subsp. japonica Proteasome subunit alpha type-7-A (Os08g0548900, LOC_Os08g43540)

Product Code	CSB-MP747735OFG
Abbreviation	Os08g0548900, LOC_Os08g43540
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YT00
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MARYDRAITV FSPDGHLFQV EYALEAVRKG NAAVGVRGTD TVVLGVEKKS TPKLQDSRSM RKIASLDTHI ALACAGLKAD ARVLINRARV ECQSHRLTVE DPVTVEYITR YIAGLQQKYT QSGGVRPFGL STLIVGFDPY TEKPALYQTD PSGTFSAWKA NATGRNSNSM REFLEKNYKD TSGKETIKLA IRALLEVVES GGKNIEIAVM THKDGLRELE EAEIDEYVAE IEAEKAAAEA AKKGAPKET
Source	Mammalian cell
Target Names	Os08g0548900
Protein Names	Recommended name: Proteasome subunit alpha type-7-A EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit D-1 20S proteasome subunit alpha-4-A
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.