



Recombinant *Oryza sativa* subsp. japonica Homeobox-leucine zipper protein HOX5 (HOX5)

Product Code	CSB-YP747762OFG
Abbreviation	HOX5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZA74
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDPGRVVFDS GVARRACPGG AQMLLFGGGG SANS GGFFRG VPAAVLGMDE SRSSSSAAGA GAKRPFTH EELLEEEYYD EQAPEKKRRL TAEQVQMLER SFEEENKLEP ERKTELARRL GMAPRQVAVW FQNRARWKT KQLEHDFDRL KAAYDALAAD HHALLSDNDR LRAQVISLTE KLQDKETSPS SATITAAQE VDQPDEHTEA ASTTG FATVD GALAAPP PGH QQPPHKDDL V SSGGTNDGD GGA AVVFDV TEGANDRLSC ESAYFADAAE AYERDCAGHY ALSSEEEDGG AVSDEGCSFD LPDAAAAAAAA MFGAAGVVHH DAADDEEAQL GSWTAWFWS
Source	Yeast
Target Names	HOX5
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX5 Alternative name(s): HD-ZIP protein HOX5 Homeodomain transcription factor HOX5 OsHox5
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.