



# Recombinant Human Zinc finger protein AEBP2 (AEBP2)

<b>Product Code</b>	CSB-BP747778HU
<b>Abbreviation</b>	AEBP2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZN18
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AAAITDMAD LEELSRLSPL PPGSPGSAAR GRAEPPPEEEEE EEEEEEEEEAE AEAVAALLLN GSGGGGGGGG GGGVGGGEAE TMSEPSESA SQAGEDEDEE EDDEEEEDDES SSSGGGEEES SAESLVGSSG GSSSDETRSL SPGAASSSSG DGDGKEGLEE PKGPRGSQGG GGGGSSSSSV VSSGGDEGYG TGGGGSSATS GGRRGSLEMS SDGEPLSRMD SEDSISSTIM DVDSTISSGR STPAMMNGQG STTSSSKNIA YNCCWDQCQA CFNSSPDLAD HIRSIHVDGQ RGGVFLVCLWK GCKVYNTPTST SQSWLQRHML THSGDKPFKC VVGGCNASFA SQGGLARHVP THFSQQNSSK VSSQPKAKEE SPSKAGMNR RKLKNKRRRS LPRPHDFFDA QTLDAIRHRA ICFNLSAHIE SLGKGHSVVF HSTVIAKRKE DSGKIKLLH WMPEDILPDV WVNESERHQL KTKVVHLSKL PKDTALLLD NIYRTMPQKR LKRTLIRKVF NLYLSKQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	AEBP2
<b>Protein Names</b>	Recommended name: Zinc finger protein AEBP2 Alternative name(s): Adipocyte enhancer-binding protein 2 Short name= AE-binding protein 2
<b>Expression Region</b>	2-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.