



Recombinant Human Protein FAM47E (FAM47E)

Product Code	CSB-EP747809HU
Abbreviation	FAM47E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZV65
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MADRRRRRLRP GTLAPVREGV NCRSRCFTKH KNGLKFPTSL HSRQLVFPRK GLDDFRKGGCP PCTGLVTQVP VEGFLPQIYH RAPQLAPKKR QIKLLKEADV LSKLSPAQQA RKAFLLEDVEA HLTPHPLALY LNLEEAMPIE LLSKVLEVLD PDRKLEDTWA YCQDTRKGMK EPTKLLKKHS TQVYLGPSKK TSVSNAGQWL YEEKPHKMDL LHENGPRPGL HENGISDIDE EFILKQFDID YETKPSHDAL HTMKLNQVPL ELKRSVGLSK LQETEFFQKL GYERKLQKPQ NPYKPKWVKM RYGAWYLNPK LWKKQRVDEP LVDPEVSHKA QEENFKKELQ EQEELLADLH GTVAFKDFIL SRGYRTPRFL ENMYIGKECK RACNKTPIKR TQA
Source	E.coli
Target Names	FAM47E
Protein Names	Recommended name: Protein FAM47E
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.