



# Recombinant Human TOM1-like protein 2 (TOM1L2)

<b>Product Code</b>	CSB-BP747812HU
<b>Abbreviation</b>	TOM1L2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZVM7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEFLLGNPFS TPVGQCLEKA TDGSLQSEDW TLNMEICDII NETEEGPKDA IRALKKRLNG NRNYREVMLA LTVLETCVKN CGHRFHILVA NRDFIDSVLV KIISPKNNPP TIVQDKVLAL IQAWADAFRS SPDLTGVVHI YEELKRKGV E FPMADLDALS PIHTPQRSVP EVDPAATMPR SQQQRTSAG SYSSPPPAPY SAPQAPALSV TGPITANSEQ IARLRSELDV VRGNTKVMSE MLTEMVPGQE DSSDLELLQE LNRTC RAMQQ RIVELISRVS NEEVTEELLH VNDDLNNVFL RYERFERYS GRSVQNASNG VLNEVTEDNL IDLGGPSPAV VSPMVGNTAP PSSLSSQLAG LDLGTESVSG TLSSLQCCNP RDGFDMFAQT RGNSLAEQRK TVTYEDPQAV GGLASALDNR KQSSEGIPVA QPSVMDDIEV WLRTDLKGGD LEEGVTSEEF DKFLEERAKA AEMVPDLSP PMEAPAPASN PSGRKKPERS EDALFAL
<b>Source</b>	Baculovirus
<b>Target Names</b>	TOM1L2
<b>Protein Names</b>	Recommended name: TOM1-like protein 2 Alternative name(s): Target of Myb-like protein 2
<b>Expression Region</b>	1-507
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.