



Recombinant Human TOM1-like protein 2 (TOM1L2)

Product Code	CSB-MP747812HU
Abbreviation	TOM1L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZVM7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEFLLGNPFS TPVGQCLEKA TDGSLQSEDW TLNMEICDII NETEEGPKDA IRALKKRLNG NRNYREVMLA LTVLETCVKN CGHRFHILVA NRDFIDSVLV KIISPKNNPP TIVQDKVLAL IQAWADAFRS SPDLTGVVHI YEELKRKGV E FPMADLDALS PIHTPQRSVP EVDPAATMPR SQQQRTSAG SYSSPPPAPY SAPQAPALSV TGPITANSEQ IARLRSELDV VRGNTKVMSE MLTEMVPGQE DSSDLELLQE LNRTC RAMQQ RIVELISRVS NEEVTEELLH VNDDLNNVFL RYERFERYS GRSVQNASNG VLNEVTEDNL IDLGGPSPAV VSPMVGNTAP PSSLSSQLAG LDLGTEVS VG TLSSLQQCNP RDGFDMFAQT RGNSLAEQRK TVTYEDPQAV GGLASALDNR QQSSEGIPVA QPSVMDDIEV WLRTDLKGGD LEEGVTSEEF DKFLEERAKA AEMVPDLSP PMEAPAPASN PSGRKKPERS EDALFAL
Source	Mammalian cell
Target Names	TOM1L2
Protein Names	Recommended name: TOM1-like protein 2 Alternative name(s): Target of Myb-like protein 2
Expression Region	1-507
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.