



Recombinant Human Inactive polyglycyclase

TTLL10 (TTLL10)

Product Code	CSB-MP747813HU
Abbreviation	TTLL10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZVT0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDHSCTRFIH RRGPPTRTRA GFKRGKRPRI QQRPRARVSG TIPASRLHPA PASQPGPCPA PGHCPVGPAH ERPMGSSQEE GLRCQPSQPD HDADGHCGPD LEGAERASAT PGPPGLLSH RPADSDDTNA AGPSAALLEG LLLGGGKPSF HSTRPGPFY IGGSNGATII SSYCKSKGWQ RIHDSRRDDY TLKWCEVKSR DSYGSFREGE QLLYQLPNNK LLTTKIGLLS TLRGRARAMS KASKVPGGVQ ARLEKDAAAP ALEDLPWTSP GYLPRQVRVLR MEEFFPETYR LDLKHREAF FTLFDETQIW ICKPTASNQG KGIFLLRNQE EVAALQAKTR SMEDDPIHHK TPFRRGPQARV VQRYIQNPPL VDGRKFDVRS YLLIACTTPY MIFFGHGYAR LTLSLYDPHS SDLGGHLTNQ FMQKKSPLYM LLKEHTVWSM EHLNRYISDT FWKARGLAKD WVFTTLKKRM QQIMAHCFLA AKPKLDCKLG YFDLIGCDFL IDDNFKVWLL EMNSNPALHT NCEVLKEVIP GVVIETLDLV LETFRKSLRG QKMLPLLSQR RFVLLHNGEA DPRPHLGGSC SLRRWPPLPT RQAKSSGPPM PHAPDQPGAR RPAPPPLVPQ RPRPPGPDLD SAHDGEPQAP GTEQSGTGNR HPAQEPSPGT AKEEREPEEN ARP
Source	Mammalian cell
Target Names	TTLL10
Protein Names	Recommended name: Inactive polyglycyclase TTLL10 Alternative name(s): Tubulin--tyrosine ligase-like protein 10
Expression Region	1-673
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.