



Recombinant *Ashbya gossypii* Vacuolar-sorting protein SNF7 (SNF7)

Product Code	CSB-EP748346DOT
Abbreviation	SNF7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q753W3
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MWSYLFSGSN AKQKKELPKK AILELREHIQ LLSKKQAHQ TQMTSQEESA RHFLSKGNKN MAKALKRKK VFESQLLKID KQVDSLEQQL YSIENANLNL ETMKAMKQGA SAMKHIHEGL DVERVDETMD EIREQVELGE EISEAISRPL YTGVEIDED ELDEELDMLA QEEVARKVME PGAAPVAAAP ATKVSLPSVP SSDMPQQDRG QPVTLGMDQE EDEDERALRE LQAEMGL
Source	E.coli
Target Names	SNF7
Protein Names	Recommended name: Vacuolar-sorting protein SNF7 Alternative name(s): Vacuolar protein-sorting-associated protein 32
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.