



Recombinant *Ashbya gossypii* Eukaryotic translation initiation factor 3 subunit I (TIF34)

Product Code	CSB-YP748380DOT
Abbreviation	TIF34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q759L2
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MRPIMLMGHE RSLTQVKYNR EGDLIFTSGK DNVASVWYAM NGERLGTLEG HNGSIWSIDV DQHTEYAVTG SADFSVKVWR VRDGSIAHSW DTRTPVRRVE FSPTGDRVLA VLDNVMNYAG AIVVFSVTRD ANNQITGFNS GLSCEILTQE GCAPVLVASW SYDGKYIVAG HQDGKISKYN GVTGECLEIK DLHKQRVSDI QFSLDRTYFL TTSRDSYANL VDVETFEVLK TYETDCPLNS GCITPLKEFV ILGGGQDARD VTTTSAREGK FEAIFYHKLK QVEIGRVDDH FGPVNYIAVS PQGTSYASGG EDGFVRLHHF DKSIFYDFKFD VEKSADAQKK VDTADR
Source	Yeast
Target Names	TIF34
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit I Short name= eIF3i Alternative name(s): Eukaryotic translation initiation factor 3 39 kDa subunit homolog Short name= eIF-3 39 kDa subunit homolog
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.