



Recombinant *Ashbya gossypii* Maintenance of telomere capping protein 2 (MTC2)

Product Code	CSB-MP748416DOT
Abbreviation	MTC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75DP8
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MVSELSLPEF EAGLSASLTC KRNFA CFATF AGEIIGALRG SRLRVQLLED FSTIGKDSVL SGELQVLVLA TMQELPPREQ NEVARWLRAL RKVHPGLVCI ATIGAHENGG PGMTNFLRRQ FWFATRGP KR ACGALRLPKR AAFTVHPAVH RYVLDVLVHI RMHRLLDHSQ AGGASSSSAE DVLDCRWLC SANHPEKTFV TPDDVQQACA WYFPMHLDVI RLPQQEASVL YGTSMKFAED LLIGLRTFLK KTGTVDNPLL LETLIVQEV L GKVVPAI
Source	Mammalian cell
Target Names	MTC2
Protein Names	Recommended name: Maintenance of telomere capping protein 2
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.