



Recombinant *Ashbya gossypii* DASH complex subunit DAM1 (DAM1)

Product Code	CSB-BP748425DOT
Abbreviation	DAM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75EQ0
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSKESVRESK AATEYRLSIS SNPGSRRSSF GGSSEHPAGS GIGNNKASNG ELREATLVDT LLLPQVQELK DSMITLDANL THMNFHESL VDLNESVSAL LYGLMCNSWC VDFPNMPHHT ARELGISKEL ARLKEEKQQL LADLQGTAQA PSLVLKEKEP NTSKQKFQLP KPPMVSTRSV VAPVRTISEE EEDDNTAASF VSNPTVMGQP PHAPPPVARS NNAKGRRRYS ILQQIRNHDL TGQKHMVANL GGTVKARAAT HPIPESAGEK RKSLAVSAVR IGNNKRLQTS RPPSGGRDVT MARKRSGTQP AILNTNSITH SSTGAVPASS AQGRRPPFR
Source	Baculovirus
Target Names	DAM1
Protein Names	Recommended name: DASH complex subunit DAM1 Alternative name(s): Outer kinetochore protein DAM1
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.