



Recombinant *Oryza sativa* subsp. japonica BURP domain-containing protein 8 (BURP8)

Product Code	CSB-EP748432OFG-B
Abbreviation	BURP8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75G46
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	GPTS PRPLWSNQGL PFRFQSKNRF SPVAAARASH AASPAELYWK IALPTSPMPG AIRDLINPAR SASQEDTMD EVSTDAVFFL EKDLFPGSKI TLHFTRGGAC AMVLLRGRAD AIPFASEKLP EILTQLSVPA GSRAAEDMRT TLAECEAALL GARDQAKHCV TSLESMVEFA AASLGTRDIR AVSTEVIGTG AAETPRQEYT VEAVKPVVSV SGGNMVTCHG MPYAYAVFGC HTTTATAYAV TLAGADGTRA EALATCHGDA FPGVAEAYER VGVAAGSVPV CHIMPLGDML WVRN
Source	E.coli
Target Names	BURP8
Protein Names	Recommended name: BURP domain-containing protein 8 Short name= OsBURP08
Expression Region	17-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.