



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 8 (TULP8)

Product Code	CSB-MP748435OFG
Abbreviation	TULP8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75HX5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSFRSIVRDV RDSFGSLR R SFEVTLAGLS GLTGHHRGKS QSTVHELCDADLIIQESRWA SLPELLRDV IRRLEASEST WPSRKDVVSC AAVCKAWREMCKEIVLSPEF CGKLTFFLSL KQPGPRDGM I QCFIKRDKSK STYHLYLCLS TAVLADSGKF LLSAKRHRKT TCTEYVISMD ADNISRSSST YIGKLRSNFLGTKFIIYDTQ PSYNGAVIPP VGRSSRRFNS KKVSPKMPSPG SYNIAQVTYE LNVLGTRGPR RMHCVMSHSIP ASSVEPGGIV PGQPEQIVPR AFEESFRSTTSFSKSSIMDR SMDFSSSRDF SSARFSDIAG GTINGDEEGQ NKERPLVLRN KAPRWHEQLQ CWCLNFRGRV TIASVKNFQL IAAPAQPPAG APTPSQPAPPEQDKIILQFG KVAKDMFTMD YRYPLSAFQA FAICLSSFDT KLACE
Source	Mammalian cell
Target Names	TULP8
Protein Names	Recommended name: Tubby-like F-box protein 8 Short name= OsTLP8 Alternative name(s): Tubby-like F-box protein 1 Short name= OsTLP1
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.