



Recombinant Chlorocebus aethiops Cytohesin-2 (CYTH2)

Product Code	CSB-EP748483DSU-B
Abbreviation	CYTH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q76MY7
Product Type	Recombinant Protein
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Purity	>85% (SDS-PAGE)
Sequence	MEDGVYEPPD LTPEERMELE NIRRRKQELL VEIQLREEL SEAMSEVEGL EANEGSKTLQ RNRKMAMGRK KFNMDPKKGI QFLVENELLQ NTPEEIARFL YKGEGLNKTA IGDYLGREE LNLAVLHAFV DLHEFTDLNL VQALRQFLWS FRLPGEAQKI DRMMEAFAQR YCLCNPGVFQ STDTCYVLSF AVIMLNTSLH NPNVRDKPGL ERFVAMNRGI NEGGDLPEEL LRNLYSIRN EPFKIPEDDG NDLTHTFNP DREGWLLKLG GRVKTWKRRW FILTDNCLYY FEYTTDKPR GIIPLENLSI REVDDPRKPN CFELYIPNNK GQLIKACKTE ADGRVVEGNH MVYRISAPTQ EEKDEWIKSI QAAVSVDPFY EMLAARKKRI SVKKKQEQP
Source	E.coli
Target Names	CYTH2
Protein Names	Recommended name: Cytohesin-2 Alternative name(s): PH, SEC7 and coiled-coil domain-containing protein 2
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.