



# Recombinant Cysteine desulfurase (sufS)

<b>Product Code</b>	CSB-MP748600EOD
<b>Abbreviation</b>	sufS
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7ADI4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTFSVDKVRA DFPVLSREVN GLPLAYLDSA ASAQKPSQVI DAEAEFYRHG YAAVHRGIHT LSAQATEKME NVRKRASLFI NARSAEELVF VRGTTEGINL VANSWGNSNV RAGDNIISQ MEHHANIVPW QMLCARVGAE LRVIPLNPDG TLQLETLPTL FDEKTRLLAI THVSNVLGTE NPLAEMITLA HQHGAKVLVD GAQAVMHPV DVQALDCDFY VFSGHKLYGP TGIGILYVKE ALLQEMPPWE GGGSMIATVS LSEGTTWTKA PWRFEAGTPN TGGIIGLGAA LEYVSAQGLN NIAEYEQNLM HYALSQLESV PDLTLYGPQN RLGVI AFNLG KHHAYDVGSF LDNYGIAVRT GHHCAMPLMA YYNVPAMCRA SLAMYNTHEE VDRLVTGLQR IHRLLG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	sufS
<b>Protein Names</b>	Recommended name: Cysteine desulfurase EC= 2.8.1.7 Alternative name(s): Selenocysteine beta-lyase Short name= SCL Selenocysteine lyase EC= 4.4.1.16 Selenocysteine reductase
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.