



Recombinant Arabidopsis thaliana Proteasome subunit beta type-7-B (PBB2)

Product Code	CSB-BP748659DOA
Abbreviation	PBB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7DLS1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	TGT TIVGLIFKDG VILGADTRAT EGPIVADKNC EKIHYPAPNI YCCGAGTAAD TEAVTDMVSS QLRLHRYQTG RDSRVVTALT LLKKHLFSYQ GHVSAALVLG GVDITGPHLH TIYPHGSTDT LPFATMGSGS LAAMSVFEAK YKEGLTRDEG IKLVAEAICS GIFNDLGSGS NVDICVITKG HKEYLRNYME PNPRTYVSSK GYSFTKKTEV LLTKITPLLE RVEIVEVAGE AMEE
Source	Baculovirus
Target Names	PBB2
Protein Names	Recommended name: Proteasome subunit beta type-7-B EC= 3.4.25.1 Alternative name(s): 20S proteasome beta subunit B-2 Proteasome component FC Proteasome subunit beta type-2
Expression Region	38-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.