



Recombinant Drosophila melanogaster Zinc finger protein GLIS2 homolog (sug)

Product Code	CSB-BP748702DLU
Abbreviation	sug
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7K0S9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDIIQKSIFN SGPHSRGIYE PPLGYFTPYN TPPYIAAYS DSGSWLADHHQ HHQQHQHQHQ QMQHIRFPT PPITPPRPIA GYGYRQRTQS VIMKARGQQD ELCRSPVEFP DDSKSCSSSS ECGTASDFVC NWTDCDRVFD TLDALAQHVT QRHAIASLTD GLYYCRWRGC QRSERGFNAR YKMLVHTRTH TKEKPHRCHL CEKSFSRAEN LKIHIRSHSG EKPYKCSFEG CQKAYSNSSD RFKHTRTHSM EKPYMCKVAG CQKRYTDPSS LRKHVKTFKH SIHLIASQPL TLPSVPCLE ASSESAFTCL PAASSVESTS SSSSARYYDD SNNEPSDYSL KPKQDAEFSP SYWLGDRQHS YLHSEDFFFVK MDVESPLDLR IHRI
Source	Baculovirus
Target Names	sug
Protein Names	Recommended name: Zinc finger protein GLIS2 homolog Alternative name(s): Protein sugarbabe
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.