



Recombinant Drosophila melanogaster Stellate protein CG33239/CG33241 (Ste:CG33239)

Product Code	CSB-YP748712DLU
Abbreviation	Ste:CG33239
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7KV15
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSSSQNNNSS WIDWFLGIKG NQFLCRVPTD YVQDTFNQMG LEYFSEILDV ILKPVIDSSS GLLYGDEKKW YGMIHARYIR SERGLIAMHR KYMRGDFGSC PNISCDRQNT LPVGLSAVWG KSTVKIHCPR CKSNFHPKSD TQLDGAMFGP SFPDIFFSML PNLTSPLDDP RT
Source	Yeast
Target Names	Ste:CG33239
Protein Names	Recommended name: Stellate protein CG33239/CG33241
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.