



Recombinant Bovine Collagen alpha-2 (IV) chain (COL4A2)

Product Code	CSB-EP749200BO-B
Abbreviation	COL4A2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SIB3
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	ISIGYLLVKH SQTdqEPMCP VGMNKLWSGY SLLYFEGQEK AHNQDLGLAG SCLARFSTMP FLYCNP GDVC YYASRNDKSY WLSTTAPLPM MPVAEEDIRP YISRCSVCEA PAVAI AVHSQ DVSIPHCPAG WRSLWIGYSF LMHTAAGDEG GGQSLVSPGS CLEDFRATPF IECNGARGTC HYYANKYSFW LTTIPEQSFQ GTPSADTLKA GLIRTHISRC QVCMKNL
Source	E.coli
Target Names	COL4A2
Protein Names	Recommended name: Collagen alpha-2(IV) chain
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes one of the six subunits of type IV collagen, the major structural component of basement membranes. The C-terminal portion of the protein, known as canstatin, is an inhibitor of angiogenesis and tumor growth. Like the other members of the type IV collagen gene family, this gene is organized in a head-to-head conformation with another type IV collagen gene so that each gene pair shares a common promoter.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.