



Recombinant Danio rerio RISC-loading complex subunit tarbp2 (tarbp2)

Product Code	CSB-EP749207DIL
Abbreviation	tarbp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SXR1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MNDEKSQDNG KRNSGYTSIE QMLAVNPGKT PISLLQEYGT RIGKTPVYDL LKAEGQAHQP NFTFRVSVGD INCTGHGPSK KAAKHKAAEA ALKMLKGGML GGIGGNGMEG DGFVGIEMEG ECPQSEMKSS SSTQQAECNP VGALQELVVQ KGWRLPEYTV TQESGPAHRK EFTMTCRVER FVEIGSGTSK KLAKRNAAAK MLSRIHDVPV DMRSSHEAEA EDDTFNMQIG GRLEGGKSKG LGCTWDSLRLN SAGEKILQLR CHPLGQSDSI DSNFCSSLRE LSEEQRFGVS YLDIEERSLS GLYQCLVELS TQPITVCHGF ASSLDAARAS AAHNALQYLK IMAGGK
Source	E.coli
Target Names	tarbp2
Protein Names	Recommended name: RISC-loading complex subunit tarbp2
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.