



Recombinant Danio rerio Zinc finger protein AEBP2 (aebp2)

Product Code	CSB-YP749208DIL
Abbreviation	aebp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SXV2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MASTTSDASE QESANPSEKN GDGSEGLRA ADPEPAVKMD EADTDEKPGG HEPTEQSPEE ATRDGGNDED GETEKEQQKS PGSADGTAER MKREQREEES VSLPAEKNSG GGGTGMKRRA SLEMMSSSSD GEPLSRMDSE DSINSTLMDT ESIPSSGRST PAMMNGLGGA SSSTKGSSYP CCWDQCHMLF NTSPDLAEHI RGVHVDGQRG GVFVCLWKGC KVYNTPSTSQ SWLQRHMLSH SGDKPFKCVV GGCNATFASQ GGLARHVPTH FSSQNSSKLS SQTKVKEESP SKAGMKNKRRK LKYKRRRSLP RPHDFFDAQT MDAIRHRAIC LNLATHIESQ GKGHSVVFHS TVIARRKEDS GKVKVLLHWT PEDILPDSWV SESDRPQLKT KVVHLSQLPQ DTAFLLDPNI YRMFF
Source	Yeast
Target Names	aebp2
Protein Names	Recommended name: Zinc finger protein AEBP2 Alternative name(s): Adipocyte enhancer-binding protein 2 homolog Short name= AE-binding protein 2 homolog
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.