



Recombinant Danio rerio RAB6-interacting golgin (gorab)

Product Code	CSB-MP749231DIL
Abbreviation	gorab
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7T320
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MAAWAGFSEE ELRKLQHNGE VANQAVSVTL GRGRKPAQTN RSRQQLQRER ALQLAAKQKE ISQSLLPDQQ LTRPPEPPAS SHPAPPAETP DEQSVPVKHE AEPVKPDPAS PVQSIPPADF KELDKQEVEL REKNRLQQQLQ WEQRIMEEKN KKRKALLTKT IAEKSKQTQA EAIKLLKIQR ELQVLDDSVS SDIGVLRKLI EQSSMDYSLA WKRFKAEAE YVAAKMDLHR KTEVKEQLTE HLCIIQQNE LRKARKLEEL MLQLELNAEE VPSPDQTQLQ DKQKEPVTEN GPAELERCT EANGSSMSKD SSITETKISQ DNQESEAADV SADGLR
Source	Mammalian cell
Target Names	gorab
Protein Names	Recommended name: RAB6-interacting golgin Alternative name(s): N-terminal kinase-like-binding protein 1 Short name= NTKL-BP1 Short name= NTKL-binding protein 1 SCY1-like 1-binding protein 1 Short name= SCYL1-BP1 Short
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.