



# Recombinant Bordetella bronchiseptica Phosphatidylserine decarboxylase proenzyme (psd)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP749823BTY   |
| <b>Abbreviation</b>      | psd   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q7WIF7  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Bordetella bronchiseptica (strain ATCC BAA-588 / NCTC 13252 / RB50)<br>(Alcaligenes bronchisepticus)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MPFKDQLFLA SQYLAPHHLV SRLMGRAADC RAPEIKNRM I ARFVRRYNVD<br>MSEALVEDPL AYASFNDFFT RALKPDARPL DDEPGAALCP ADGAISQIGA<br>IDNGRIFQAK GHSFGLTDLL GGDAERAAPF AGGQFATIYL SPRDYHRVHM<br>PLAGTLREMV HVPGR LFSVN PLTARSVPEL FARNERVA CL FDTEHGPMAL<br>VLVGAMIVAS IETVWAGLVT PHKRQVRSVR YDAAARAPIH LDKGAEMGRF<br>KLG  |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | psd   |
| <b>Protein Names</b>     | Recommended name: Phosphatidylserine decarboxylase proenzyme EC=<br>4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine<br>decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain   |
| <b>Expression Region</b> | 1-253   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at<br>4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the<br>contents to the bottom. Please reconstitute protein in deionized sterile water to a<br>concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final<br>concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final<br>concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,<br>storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.