



Recombinant Staphylococcus aureus Pyruvate dehydrogenase E1 component subunit beta (pdhB)

Product Code	CSB-MP750025SKW
Abbreviation	pdhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GAC0
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain MSSA476)
Purity	≥85% (SDS-PAGE)
Sequence	MAQMTMVQAI NDALKTELKN DQDVLIFGED VGVNGGVFRV TEGLQKEFGE DRVFDTPLAE SGIGGLAMGL AVEGFRPVME VQFLGFVFEV FDAIAGQIAR TRFRSGGTKT APVTIRSPFG GGVHTPELHA DNLEGILAQS PGLKVVIPSG PYDAKGLLIS SIRSNDPVVY LEHMKLYRSF REEVPEEEYT IDIGKANVKK EGNDISIITY GAMVQESMKA AEELEKDGYS VEVIDLRTVQ PIDVDTIVAS VEKTGRAVVV QEAQRQAGVG AAVVAELSER AILSLEAPIG RVAAADTIYP FTQAENVWLP NKNDIIEKAK ETLEF
Source	Mammalian cell
Target Names	pdhB
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit beta EC= 1.2.4.1
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.