



Recombinant *Xenopus laevis* Cysteine protease ATG4A (atg4a)

Product Code	CSB-BP750206XBE
Abbreviation	atg4a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GPU1
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MSDPVS DYLY KYENEPEYLD LEELPDSDEP VYILGKQYDT KTDKCDLQSD IVSRLWFTYR KKFSPIGGTG PSSDTGWGCM LRCGQMMLAQ ALVCQHLGRD WRWEKHKNHP EYQQILQCF LDRKDCCYSI HQMAQMGVGE GKSIGEWFGP NTVAQVLKKL ALFDEWNSLA VYVSM DNTVV VEDIKTMCKY QPQSCSMAQA ASHQSTWSRC RDTSGHCSGW RPLLLVPLR LGINHINPVY VDAFKACFKM PQSLGALGGK PNHAYYFIGF SGDEIYLDP HTTQTFVDTE EAGTVQDQTY HCQKGPNSMK VLNLDP SVAL GFFCKDENDF NNWCEVIEKE ILKHQSLRMF ELTPKHPPHW PPFIPPTKPE VTTTGAE LIE STDKLF DVEE EFEILSV
Source	Baculovirus
Target Names	atg4a
Protein Names	Recommended name: Cysteine protease ATG4A EC= 3.4.22.- Alternative name(s): Autophagy-related protein 4 homolog A
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.