



Recombinant Mouse Dual specificity protein phosphatase CDC14A (Cdc14a)

Product Code	CSB-BP750214MO
Abbreviation	Cdc14a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GQT0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAAESGELIG ACEFMKDRLY FATLRNRPKS TINIHYSID EELVYENFYA DFGPLNLAMV YRYCCKLNKK LKSYLSRKK IVHYTSFDQR KRANAFLIG AYAVIYLKKT PEEAYRALLS GSNPPYLPFR DASFGNCTYN LTVLDCLQGI RKGLQHGGFFD FETFDAEEYE HYERVENGDF NWIVPGKFLA FSGPHPKSKI ENGYPLHAPE AYFPYFKKNN VTTIVRLNKK IYEAKRFTDA GFEHYDLFFI DGSTPSDNIV RRFLNICENT EGAIHVHCKA GLGRTGTLIA CYVMKHYRFT HAEIIAWIRI CRPGSIIGPQ QHFLKEKQAS LWVQGDIFRS KLKNRPSSEG SITKIISTLD DMSIGANLSK LQSTERIGEN NFEDEDMEIK NNVTQGDKLR ALKSQRHPRS SPSCAFRSDD MKGHQRAMAQ TFRLLSSSPQP TMSTMKTSKV CLSPSVTAKK ISRGLSSGA NIRSFSINSR LASSLGNLNA GTEEPETKKT TSLTKAAFIA SPFTSFLNGS TQTPGRNYPE LNNNQYTRSS NSNSSSSSSG LGGNLNSSPV PQSAKPEEHT TILRPSFPGS LSSSSVRFLS RSIPSLQSEY VHY
Source	Baculovirus
Target Names	Cdc14a
Protein Names	Recommended name: Dual specificity protein phosphatase CDC14A EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): CDC14 cell division cycle 14 homolog A
Expression Region	1-603
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the dual specificity protein tyrosine phosphatase family. It is highly similar to Saccharomyces cerevisiae Cdc14, a protein tyrosine phosphatase involved in the exit of cell mitosis and initiation of DNA replication, suggesting a role in cell cycle control. This protein has been shown to interact



with, and dephosphorylate tumor suppressor protein p53, and is thought to regulate the function of p53. Alternative splicing of this gene results in several transcript variants encoding distinct isoforms.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.