



Recombinant *Oryza sativa* subsp. japonica Cysteine proteinase inhibitor 3 (Os05g0409300, LOC_Os05g33880)

Product Code	CSB-BP750312OFG
Abbreviation	Os05g0409300, LOC_Os05g33880
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6I570
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	FHLGG DESVLVRGML AAIRREQAEA EDAARFAVAE YNKNQGAELE FARIVKAKRQ VVTGTLHDLM LEVVDSGKKS LYSAKVWVKP WLDFKAVVEF RHVGDSQSQS ATAADDNAGQ DTADPTVASR NDLHNTENNK VSVVLSTFSQ TYSV
Source	Baculovirus
Target Names	Os05g0409300
Protein Names	Recommended name: Cysteine proteinase inhibitor 3 Alternative name(s): Oryzacystatin III Short name= OC-III Oryzacystatin-3
Expression Region	36-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.