



Recombinant *Oryza sativa* subsp. *japonica* BURP domain-containing protein 2 (BURP2)

Product Code	CSB-EP750313OFG
Abbreviation	BURP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6I5W0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ASPAEMYWK IALPTSPMPG AIRDLISPAS SAASASKDKE DTVGSVFFLE KDLFPGSKMT LHFTRATAGA ALLPRGRADS VPFASEKLPE ILSQLSIPAG SPTADAMRST LAVCEAARIA SETAPKHKHY CATSLESMVE LVASSLGTRD VHAVSTEVVN RAGPTPRQAY RVEAVRPVPV PGGDMVACHR MPYAYAVFGV HGIKGAAYTV TLAGADGTMA EAVAACHGDV DGHGVAVAEA YKRLGVAPGK VAVCHFLPQD DMLWVRN
Source	E.coli
Target Names	BURP2
Protein Names	Recommended name: BURP domain-containing protein 2 Short name= OsBURP02
Expression Region	22-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.