



Recombinant Mouse Tropomyosin alpha-4 chain (Tpm4)

Product Code	CSB-BP750349MO
Abbreviation	Tpm4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IRU2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AGLNSLEAV KRKIQALQQQ ADDAEDRAQG LQRELDGERE RREKAEGDAA ALNRRIQLE EELDRAQEQL ATALQNLEEA EKAADESERG MKVIENRAMK DEEKMEILEM QLKEAKHITD EADRYEEVA RKLVILEGEL KRAEERAIEVS ELKCGDLEEE LKNVTNNLKS LEAASEKYSE KEDKYEEEEIK LLSDKLKEAE TRAEFAERTV SKLEKTIDDL EEKLAQAKEE NVGLHQTLTDQ TLNELNCI
Source	Baculovirus
Target Names	Tpm4
Protein Names	Recommended name: Tropomyosin alpha-4 chain Alternative name(s): Tropomyosin-4
Expression Region	2-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the tropomyosin family of actin-binding proteins involved in the contractile system of striated and smooth muscles and the cytoskeleton of non-muscle cells. Tropomyosins are dimers of coiled-coil proteins that polymerize end-to-end along the major groove in most actin filaments. They provide stability to the filaments and regulate access of other actin-binding proteins. In muscle cells, they regulate muscle contraction by controlling the binding of myosin heads to the actin filament. Multiple transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.