



Recombinant Human Grainyhead-like protein 2 homolog (GRHL2)

Product Code	CSB-EP750350HU-B
Abbreviation	GRHL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ISB3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSQESDNNKR LVALVPMPD PPFNTRRAYT SEDEAWKSYL ENPLTAATKA MMSINGDEDS AAALGLLYDY YKVPDRKRL SVSKASDSQE DQEKRNCLGT SEAQSNLSGG ENRVQLKTV PVNLSLNQDH LENSREQYS ISFPESSAI PVSGITVKA EDFTPVFMAP PVHYPRGDGE EQRVVIFEQT QYDVPSLATH SAYLKDDQRS TPDSTYSESF KDAATEKFRS ASVGAEEMY DQTSSGTFQY TLEATKSLRQ KQGEGPMTYL NKGQFYAITL SETGDNKCFR HPISKVRSV MVFSEDKNR DEQLKYWKYW HSRQHTAKQR VLDIADYKES FNTIGNIEE AYNAVSFTWD VNEEAKIFIT VNCLSTDFSS QKGVKGLPLM IQIDTYSYNN RSNKPIHRAY CQIKVFCDKG AERKIRDEER KQNRKKGKGQ ASQTQCNSSS DGKLAAPLQ KKS DITYFKT MPDLHSQPVL FIPDVHFANL QRTGQVYYNT DDEREGGSVL VKRMFRPMEE EFGPVPKQM KEEGTKRVLL YVRKETDDVF DALMLKSPTV KGLMEASEK YGLPVEKIAK LYKKS KKGIL VNMDDNIIEH YSNEDTFILN MESMVEGFKV TLMEI
Source	E.coli
Target Names	GRHL2
Protein Names	Recommended name: Grainyhead-like protein 2 homolog Alternative name(s): Brother of mammalian grainyhead Transcription factor CP2-like 3
Expression Region	1-625
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.