



# Recombinant Human Uncharacterized protein KIAA1841 (KIAA1841)

<b>Product Code</b>	CSB-MP750733HU
<b>Abbreviation</b>	KIAA1841
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NSI8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRGYSENNN FLNNNNQMV L DMILYPLIGI PQTINWETIA RLVPG LTPKE CAKRFDELKS SGSSPVDNQY NSLMAAGESP VETLATYIKS SLLDIHGFEQ ETPVGHDAVS KTGRHSIAST RNCSSSEENC TTHNGGEMTE ESEGNMVIH VCDEAKNLKE DFTCPDLLI SEMKYFAEYL SMDAQRWEEV DISVHCDVHI FNWLIKYIKR NTKENKDC EM PTLEPGNVIS ILISSEFLKM DSLVEQCIQY CHKNMNAIVA TPCNMNCINA NLLTRIADLF SHNEVDDLKD KKDKFKSKLF CKKIERL FDP EYLNPD SRSN AATLYRCCLC KLLTKETER RIPCIPGKIN VDRRG NIVYI HIRDKTWDVH EYLN SLFEEL KSWRDVYWRL WGTINWLTCS RCYQAFLCIE FSHCQYHSET VVYPTAASSL NTVGTGIYPC CNQKVLRFDP TQLTKGCKVR DHMVTLRDQG EGGDLPSCPT ARMLDDLHKY RDVIVVPFSK DTVSDVGVGL CDEK GIECDV LLEPNTPWGP KTGELNAFLS LKNWTLQLKQ QSLFSEEEY T TGSEVTEDE VGDEEEVSKK QRKKEKPKKF TRQPKKQVSS PCAQRKEKAL EKSASRDVSP FVMSMQKNKW DATRSLRFNQ DAQREDDQRR MTEITGHLIK MRLGDLDRVK SKEAKEFAGG IYSRLEAIK ASVPVSARQS SSEKNTRSKS RFGQGRPA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KIAA1841
<b>Protein Names</b>	Recommended name: Uncharacterized protein KIAA1841
<b>Expression Region</b>	1-718
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.