



Recombinant *Xenopus tropicalis* Transcription factor Sox-2 (sox2)

Product Code	CSB-MP750756XBF
Abbreviation	sox2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NVN0
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MYNMMETDLK PPAPQQASGG NSNSGSNNQS KNSPDRVKRP MNAFMVWSRG QRRKMAQENP KMHNSEISKR LGAEWKLLSE AEKRPFIDEA KRLRALHMKE HPDYKYRPRR KTKTLMKKDK YTLPGGLLAP GANPMTSGVG ASLGAGVNQR MDTYAHMNGW TNGGYGMMQE QLGYPQHPGL SAHNAPQMMP MHRYDVSALQ YNSMSSSQTY MNGSPTYSMS YSQQGAPGMS LGSMGSVVKS ESSSSPPVVT SSSHSRIPCQ AGDLRDMISM YLPGAEPPEP AAQSRLHMSQ HYQSASVAGT AINGTLPLSH M
Source	Mammalian cell
Target Names	sox2
Protein Names	Recommended name: Transcription factor Sox-2 Alternative name(s): SRY (sex determining region Y)-box 2
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.