



Recombinant Danio rerio Cysteine dioxygenase type 1 (cdo1)

Product Code	CSB-MP750762DIL
Abbreviation	cdo1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NWZ9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEQTEVMKPE TLEDLIKTLH QIFQSDSINV EEVQNLMSY QSNPQDWMKF AKFDQYRYTR NLVDEGNGKF NLMILCWGEG HGSSIHDHTD SHCFLKLLQG QLKETLFDWP DRKLQSGMKP RGQSVLQENQ CAYINDSLGL HRVENVSHT PAVSLHLYSP PFQSCRTFDQ RTGHHNTVKM TFWISKYGERT PYELSVSQEN N
Source	Mammalian cell
Target Names	cdo1
Protein Names	Recommended name: Cysteine dioxygenase type 1 EC= 1.13.11.20 Alternative name(s): Cysteine dioxygenase type I Short name= CDO Short name= CDO-I
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.