



Recombinant Human Protein FAM71F2 (FAM71F2)

Product Code	CSB-MP750765HU
Abbreviation	FAM71F2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NXP2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSKIRGLPPE VREPGPGVEL GVENGLLCQL IHSPEFNLFS NSVVFESNFI QTHVPEADFDQ VTKPGNWRDV CEGSATVILG VTSSVPSLPL PNVLLMANVT WPQGPFTTWS TPGDAPVINL SROLLPKYVE LRIYDRLQRI LRVRTVTEKI YYLKLHEKHP EIVFQFWVRL VKILQKGLSI TTKDPRIKFT HCLVPKMPTN STETTPENSL LSSPQPSEPL VLLAAEQTSG SFSQLSGKPQ LTADRNNDTA IEIDNCSSYK IPSPVASPIN LNIPMRAALS HSLWEQEDWN EHLLQVHIAS YLGEHFLGA
Source	Mammalian cell
Target Names	FAM71F2
Protein Names	Recommended name: Protein FAM71F2 Alternative name(s): Protein FAM137B
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.