



# Recombinant Mouse Cell cycle checkpoint protein RAD17 (Rad17)

<b>Product Code</b>	CSB-YP750766MO
<b>Abbreviation</b>	Rad17
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NXW6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSETFLRPKV SSTKVTDWVA PAFDDFEANT AITTITASSL TFSNSSHRRK YLPSTLESNR LSARKRGRLS LEQTHGLETS RERLSDNEPW VDKYKPETQH ELAVHKKKIE EVETWLKAQV LEVKPKQGGG VLLITGPPGC GKTTTTIKILS KELGIQVQEW VNPILPDFQK DDYKELLSLE SNFSVVPYQS QIAVFNDFLL RATKYSKLQM LGDDLTTDDK IILVEELPNQ FYRDPNALHE ILRKHVQIGR CPLVFIVSDS VSGDNNQRLL FPRNIQEECS VSNISFNPVA PTIMMKFLNR IVTIEASKNG EKIIVPNKTS LELLCQGCSG DIRSAINSLQ FSSSKGENSS WSKKKRMSLK SDAAISKSKQ KKKHNSTLEN QEIQAIKKD VSLFLFRALG KILYCKRAPL TELDSPRLPA HLSEHDRDTL LVQPEEIVEM SHMPGDFNL YLHQNYIDFF AEVDDLVPAS EFLSFADILG GDWNTRSLLR EYSTSVATRG VMHSNKARGF AHCQGGSSFR PLHKPQWFLI QKKYRENCLA AKALFVDFCL PALCLQTQLL PYLALLTIPM RNKAQISFIQ DVGRLPLKRS FGRLKMEALT DRELGLIDPD SGDESPHSGG QPAQEAPGEP AQAAQNADPE TWSLPLSQNS GSDLPASQPQ PFSSKVDMEE EEEEEEDIII EDYDSEET
<b>Source</b>	Yeast
<b>Target Names</b>	Rad17
<b>Protein Names</b>	Recommended name: Cell cycle checkpoint protein RAD17
<b>Expression Region</b>	1-688
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is highly similar to the gene product of Schizosaccharomyces pombe rad17, a cell cycle checkpoint gene required for cell cycle arrest and DNA damage repair in response to DNA damage. This protein shares strong similarity with DNA replication factor C (RFC), and can form a complex with RFCs. This protein binds to chromatin prior to DNA damage and is



phosphorylated by the checkpoint kinase ATR following damage. This protein recruits the RAD1-RAD9-HUS1 checkpoint protein complex onto chromatin after DNA damage, which may be required for its phosphorylation. The phosphorylation of this protein is required for the DNA-damage-induced cell cycle G2 arrest, and is thought to be a critical early event during checkpoint signaling in DNA-damaged cells. Eight alternatively spliced transcript variants of this gene, which encode four distinct proteins, have been reported. Two pseudogenes, located on chromosomes 7 and 13, have been identified.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.