



Recombinant Danio rerio Pre-mRNA-splicing factor RBM22 (rbm22)

Product Code	CSB-YP750778DIL
Abbreviation	rbm22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NZZ9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MATSLGSNTY NRQNWEDSDF PILCQTCLGE NPYIRMTKEK FGKECKICAR PFTVFRWCPG VRMRFKKTEV CQTCSKMKNV CQTCLLDLEY GLPIQVRDTG LSVKDEVPRS DVNKEYYTQN MEREIANS DG TRPVGLLGKA PSSSDMLLKL ARTTPYYKRN RPHICSFVVK GECKRGEECP YRHEKPTDPD DPLADQNIKD RYYGINDPVA NKLLMRASTM PRLDVPDDKS ITTLYIGGLG ENVTDSELRN HFYQFGEIRT ITIVQRQQCA FIQFATRQAA ETAAEKSFNK LIINGRRLNV KWGRSQAARG KGEKDGVTES GIRLEPVPGL PGALPPPPVS DEDASTNYFN LAPTPSPAVM NLGLPPPPGV TMPPPPGFGP PMFHTMDPMA PPVPPPMALR PPGQVHYPSSQ DPQRMGAHAS RHGGP
Source	Yeast
Target Names	rbm22
Protein Names	Recommended name: Pre-mRNA-splicing factor RBM22 Alternative name(s): RNA-binding motif protein 22
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.