



Recombinant Human Parafibromin (CDC73)

Product Code	CSB-MP750783HU
Abbreviation	CDC73
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P1J9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ADVLSVLRQ YNIQKKEIVV KGDEVIFGEF SWPKNVKNTNY VVWGTGKEGQ PREYYTLDSI LFLNNVHLS HPVYVRRRAAT ENIPVRRPD RKDLLGYLNG EASTSASIDR SAPLEIGLQR STQVKRAADE VLAEAKKPRI EDEECVRLDK ERLAARLEGH KEGIVQTEQI RSLSEAMSVE KIAAIKAKIM AKKRSTIKTD LDDITALKQ RSFVDAEVDV TRDIVSRERV WRTRTTILQS TGKNFSKNIF AILQSVKARE EGRAPEQRPA PNAAPVDPTL RTKQPIPAAY NRYDQERFKG KEETEGFKID TMGTYHGMLT KSVTEGASAR KTQTPAAQPV PRPVSQARPP PNQKKSRTTP IIIIPAATTS LITMLNAKDL LQDLKFVPSD EKKKQGCQRE NETLIQRRKD QMQPGGTAIS VTPYRVVDQ PLKLMQDWD RVVAVFVQGP AWQFKGWPWL LPDGSPVDIF AKIKAFHLKY DEVRLDPNVQ KWDVTVLELS YHKRHLDPRV FLRFWETLDR YMVKHKSHLR F
Source	Mammalian cell
Target Names	CDC73
Protein Names	Recommended name: Parafibromin Alternative name(s): Cell division cycle protein 73 homolog Hyperparathyroidism 2 protein
Expression Region	2-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a tumor suppressor that is involved in transcriptional and post-transcriptional control pathways. The protein is a component of the the PAF protein complex, which associates with the RNA polymerase II subunit POLR2A and with a histone methyltransferase complex. This protein appears to facilitate the association of 3 mRNA processing factors with actively-transcribed chromatin. Mutations in this gene have been linked to hyperparathyroidism-jaw tumor syndrome, familial isolated hyperparathyroidism, and parathyroid carcinoma.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.