



Recombinant Danio rerio Serine/threonine-protein phosphatase 4 regulatory subunit 2-A (ppp4r2a)

Product Code	CSB-MP750830DIL
Abbreviation	ppp4r2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P964
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEIDTLLGAF RDFEKKGEKE SCPILDQFLS HVAKTGETMI QWSQFKSYFL FKLEKVMDDF KASCPEQRGP ANPNVEYIPF EEMKQRILKI VNGYNGIPFT IQRLCELLTE PKRNYAGTDK FLRGVEKNVM VVSCVYPTSE KKGSSCVNRM NGVMFPGNTS AFPDRNVNGP GTPRPLNRPK HSLSSNVATN GLPDSTESKE QASEQSERTV NESSASEAES HSGAVKSKHR DDEDATHAET HEAKRLKFDK EEEEEEEDDEE EDEDGDEIKK ELDEPHSPCT SVAESSSDVP QSSTDVTAEV KDEDQEPSST QSEVVENGVD RSTSESDPDP SHKATGSESD PKEQQAEEEE EEESAEAQET EETNDPVSSS SNNSSDEGVS SAETPSASPS SSTEPAEGS VTAEITSDNS ETADDNMEQD
Source	Mammalian cell
Target Names	ppp4r2a
Protein Names	Recommended name: Serine/threonine-protein phosphatase 4 regulatory subunit 2-A
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.