



Recombinant Mouse Ubiquitin carboxyl-terminal hydrolase MINDY-2 (Mindy2)

Product Code	CSB-MP750855MO
Abbreviation	Fam63b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PDI6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MENSPDSPQP LELGVAAGRV SPPEGRRRGG REAEDGPAGR AVDSGGQGAA AAAARSSLGD PTSPSQLGCG AGSDLKDGAS SSPAASEVPS RGQHKVTASP ELAEAAAGRG SGPVGDGTGTC RVEQAAEEPS STGAPSSSCS EPSPPGDSPS LDSLESFSNL HSFSSSEFN SEEGAETRVV EDVEEGAAGP PRAAPLCKEE EEDPAQVLA SKERFPQSV YHIKWIQWKE ENTPIITQNE NGPCPLLAIL NVLLLAWKVK LPPMMEITA EQLMEYLG DY MLEAKPKEIS EIQRVNYEQN MSDAMAILHK LQTGLDVNVR FTGVRVFEYT PECIVFDLLD IPLYHGWLVD PQIDDIVKAV GNCSYNQLVE KIISCKQSDN SQLVSEGFVA EQFLNNTATQ LTYHGLCELT STVQEGELCV FFRNNHFSTM TKYKGQLYLL VTDQGFLTEE KIVWESLHNV DGDGNFCDSE FHLRPPSDPE TVYKGGQQDQI DQDYLMALSL QQEQQSQEIN WEQIPEGISD LELAKKLQEE EDRRASQYYQ EQEQAQAVVT TTPSTQAQQ GQPAQASPSS IKQPGNSERK RKEPREKDKE KEKEKNSCVI L
Source	Mammalian cell
Target Names	Mindy2
Protein Names	Recommended name: Protein FAM63B
Expression Region	1-601
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.