



Recombinant Danio rerio Neuroguidin (ngdn)

Product Code	CSB-MP750864DIL
Abbreviation	ngdn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PFJ1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAATAVGNDV IDRDLPTAVR LLNSLTEQIV SVTSHVRELI KKVREKAYQT SKGLSFLDLR YHLLLFYQLD ITHLISLKTE GESLKDNSAI HRLVTIRTVL EKMRPLDQKL KYQIDKLVRT AVTGSLSEND PLHFRPNPQS LVSKLSESED SDDDGVGKKT KEQKEPSGGR RYVPPRIAPM HYDGDLTEAD RQKERVEKQK RAALRSSVIQ ELRQQYSDAP EEIRDRRDFQ TDRQGREELN RKNYEESMMV RLSTTRDQKL RKRGMGMGMTS QLSNITRFSD ITALTGGEVQ DIGNPKPKKK KIIKKGKKKV FRKRK
Source	Mammalian cell
Target Names	ngdn
Protein Names	Recommended name: Neuroguidin Alternative name(s): EIF4E-binding protein
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.