



Recombinant Human Protein aurora borealis (BORA)

Product Code	CSB-EP750873HU-B
Abbreviation	BORA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PGQ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGDVKESKMQ ITPETPGRIP VLNPFESPSD YSNLHEQTLA SPSVFKSTKL PTPGKFRWSI DQLAVINPVE IDPEDIHRQA LYLSHSRIDK DVEDKRQKAI EEFFTKDVIV PSPWTDHEGK QLSQCHSSKC TNINSDSPVG KKLTIHSEKS DAACQTLLSL PVDFNLENIL GDYFRADEFA DQSPGNLSSS SLRRKLFLDG NGSISDSLPS ASPGSPHSGV QTSLEMFYSI DLSPVKCRSP LQTPSSGQFS SSPIQASAKK YSLGSITSPS PISSPTFSPI EFQIGETPLS EQRKFTVHSP DASSGTNSNG ITNPCIRSPY IDGCSPIKNW SPMRLQMYSG GTQYRTSVIQ IPFTLETQGE DEEDKENIPS TDVSSPAMDA AGIHLRQFSN EASTHGTHLV VTAMSVTQNNQ SSASEKELAL LQDVEREKDN NTVDMVDPIE IADETTWIKE PVDNGSLPMT DFVSGIAFSI ENSHMCMSPS AESSVIPCES SNIQMDSGYN TQNCGSNIMD TVGAESYCKE SDAQTCEVES KSQAFNMKQD HTTQRCWMKT ASPFQCSSP
Source	E.coli
Target Names	BORA
Protein Names	Recommended name: Protein aurora borealis Short name= HsBora
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.