



# Recombinant Human Protein aurora borealis (BORA)

<b>Product Code</b>	CSB-YP750873HU
<b>Abbreviation</b>	BORA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6PGQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGDVKESKMQ ITPETPGRIP VLNPFESPSD YSNLHEQTLA SPSVFKSTKL PTPGKFRWSI DQLAVINPVE IDPEDIHRQA LYLSHSRIDK DVEDKRQKAI EEFFTKDVIV PSPWTDHEGK QLSQCHSSKC TNINSDSPVG KKLTIHSEKS DAACQTLLSL PVDFNLENIL GDYFRADEFA DQSPGNLSSS SLRRKLFLDG NGSISDSLPS ASPGSPHSGV QTSLEMFYSI DLSPVKCRSP LQTPSSGQFS SSPIQASAKK YSLGSITSPS PISSPTFSPI EFQIGETPLS EQRKFTVHSP DASSGTNSNG ITNPCIRSPY IDGCSPIKNW SPMRLQMYSG GTQYRTSVIQ IPFTLETQGE DEEDKENIPS TDVSSPAMDA AGIHLRQFSN EASTHGTHLV VTAMSVTQNNQ SSASEKELAL LQDVEREKDN NTVDMVDPIE IADETTWIKE PVDNGSLPMT DFVSGIAFSI ENSHMCMSPS AESSVIPCES SNIQMDSGYN TQNCGSNIMD TVGAESYCKE SDAQTCEVES KSQAFNMKQD HTTQRCWMKT ASPFQCSSP
<b>Source</b>	Yeast
<b>Target Names</b>	BORA
<b>Protein Names</b>	Recommended name: Protein aurora borealis Short name= HsBora
<b>Expression Region</b>	1-559
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.