



Recombinant Danio rerio Pachytene checkpoint protein 2 homolog (trip13)

Product Code	CSB-EP750876DIL
Abbreviation	trip13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PH52
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDISDHHMTD VGIIRPDVHI EVHVKSQSTA KRSDVIRTHVL SLLDRHSTVL RSFKWMDFDN EFLTKNVESV TIADVTGPKL VDLKVHNLICI HIFTLNDDSP STLNLEEEEE LSAANLWLLP AVEFHGVWES LIYEEGIKTQ LLDYVSTTIF FSDKNVDSNL IAWNRVLLH GPPGTGKTSL CKGLAQKLSI RLSDRYAHSQ FVEINSHSLF SKWFSESGKL VTKMFQKIQE LIDDKDALVF VLIDEVESLT AARSAAQAGT EPSDAIRVVN SVLTQLDQIK RHPNVVILTT SNVTEKIDLA FVDRADIKQY IGPPSAKAIF NIYLSLEEL MKRQIIYPRQ QLVSLEELET MFMESDVTR LSLCLMNISQ RSVGLSGRAL RKIPFLAHAL YGKTSTMTLK GFLTAMEHAV NKQRREQASL VNCV
Source	E.coli
Target Names	trip13
Protein Names	Recommended name: Pachytene checkpoint protein 2 homolog Alternative name(s): Thyroid hormone receptor interactor 13 homolog Thyroid receptor-interacting protein 13 homolog Short name= TR-interacting protein 13 homolog Short name=
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.