



# Recombinant Human Tripartite motif-containing protein 65 (TRIM65)

<b>Product Code</b>	CSB-BP750884HU
<b>Abbreviation</b>	TRIM65
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6PJ69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AAQLLEEKL TCAICLGLYQ DPVTLPCGHN FCGACIRDWW DRCGKACPEC REPFDPGAEL RRNVALSGVL EVVRAGPARD PGPDPGPGPD PAARCPRHGR PLELFCRTEG RCVCSVCTVR ECRLHERALL DAERLKREAQ LRASLEVTQQ QATQAEGQLL ELRKQSSQIQ NSACILASWV SGKFSSLLQA LEIQHTTALR SIEVAKTQAL AQARDEEQRL RVHLEAVARH GCRIRELLEQ VDEQTFLQES QLLQPPGPLG PLTPLQWDED QQLGDLKQLL SRLCGLLLEE GSHPGAPAKP VDLAPVEAPG PLAPVPSTVC PLRRKLWQNY RNLTDFPVSA NRHFYLSRQD QQVKHCRQSR GPGGPGSFEL WQVQCAQSFQ AGHHYWEVRA SDHSVTLGVS YPQLPRCRLG PHTDNIGRGP CSWGLCVQED SLQAWHNGEA QRLPGVSGRL LGMDLDSLASG CLTFYSLEPQ TQPLYTFHAL FNQPLTPVFW LLEGRTLTLIC HQPGAVFPLG PQEEVLS
<b>Source</b>	Baculovirus
<b>Target Names</b>	TRIM65
<b>Protein Names</b>	Recommended name: Tripartite motif-containing protein 65
<b>Expression Region</b>	2-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.